

# PHOTOFACIALS

*The Palomar Rejuvelux Process  
clears sun damage and visible veins for a Rejuvenated Appearance*

## *The RejuveLux Process*

The RejuveLux Process uses pulses of concentrated light to clear:

- \* **Sun Damage:** Wide areas of pigmentation on the face, chest, and arms.
- \* **Sun Spots:** Individual pigmented lesions, also known as age spots.
- \* **Rosacea:** General inflammation of the cheeks and face.
- \* **Spider Veins:** Individual vessels, which can be red, blue, or purple in color.
- \* **Hair Removal:** Rapid hair reduction with less discomfort and less treatment time.

A hand piece will be placed gently over your skin while pulsing. Each pulse will send a beam of light through your skin. A flashing light, an audible tone, and a mild sensation on the skin.

## *Pulse Light versus Laser Treatment*

Pulsed light systems and lasers are very similar, but Palomar Pulsed Light Systems use a much larger treatment window than most lasers. This spreads light over a wider area of skin and allows faster coverage of the skin and quicker treatments as a result.

It can treat any area of the body, except directly around the eyes. The face, neck, chest, arms, back and legs can all be treated.

## *Treatment of Choice*

This office has chosen a Palomar Pulsed Light System because it offers the most versatile, effective, and comfortable aesthetic treatments available.



### **Technical Specifications**

Pulsewidth	10/20/40/60/100 ms (Selectable)
Wavelength	Dependent on handpiece
Speed	Up to 7.4 cm <sup>2</sup> /sec (Coverage Rate), 1 Hz
Electrical	(Pulse Rate)
Req.	100-240 V, 50/60 Hz
Dimensions	23.5" (60 cm) L x 11.9" (30 cm) H x 17.8"
Weight	(45 cm) D 63 lbs. (29 kg)



**(800) 660-6162**  
[www.primedical.net](http://www.primedical.net)