

# Hoya Con Bio Medlite C6

## Faster Red, Green, Blue, and Black Tattoo Removal Phenomenal Results on Pigmented Lesions Non-Ablative Treatment for Acne Scars

The Medlite C6 Q-switched Nd:YAG does an exceptional job of clearing blue, black, and red ink with the 1064 and frequency-doubled 532 nm wavelengths. Green and sky blue inks are treated with the MultiLite dye-impregnated handpieces (585 nm and 650 nm). With spot sizes up to 8 mm, tattoos are treated quickly and effectively, with virtually no pinpoint bleeding. The beam profile lends itself to less dermal damage than with other Nd:YAG lasers, resulting in a lower incidence of scarring.

Epidermal and dermal pigmented lesions can be treated quickly and effectively, with minimal or no incidence of hypopigmentation. The 1064 wavelength also treats wrinkles and acne scars via non-ablative resurfacing.

Treating vascular lesions is easy! Physicians treating vascular lesions can now choose the IPL, Candela VBeam or the KTP Aura to treat their patients' spider veins, port wine stains, rosacea and telangiectasia.



For more information, please contact:

Rel Rodriguez  
Territory Manager  
Southern California - Aesthetic Division  
949.395.1485

Rick Hallagan  
Regional Sales Manager  
California  
916.952.7989

### **PRI Medical Technologies, Inc.**

10939 Pendleton Street  
Sun Valley, CA 91352  
www.primedical.net

**800.660.6162**

